

Lab 24, Waves, Ripple Tank Action

Name _____ Period _____

Purpose: To investigate the properties of waves in water.

1. Fill tank to the plastic frame with tap water (1cm deep).
2. Use a pencil for a point source wave generator.
3. Use a metal bar for a wave front generator.
4. Sketch the pattern of waves observed from a
 - a) point source.
 - b) wave front.
5. Sketch the pattern of wave interaction from a double point source.
6. Repeat the above with different frequencies.
7. Repeat the above with different amplitudes.

8. Reflection:

Sketch the patterns of a wave front reflected from

- a) plane reflector
- b) concave reflector.
- c) convex reflector.

9. Refraction:

Sketch the wave interactions with a shallow place (glass plate) with

- a) a direct hit.
- b) an angle attack.

10. Repeat with a different

- a) frequency
- b) amplitude.

11. Diffraction:

Sketch the wave action with a

- a) single slit
- b) double slit

12. Repeat with different

- a) frequency
- b) amplitude

13. Doppler Effect by moving

- a) the source
- b) the reflector.

14. The Grand Critique of wave interactions.