

Lab 26, Photometry

Name _____ Period _____

**** Please Bring two different candles for this Lab ****


Set up the Photometry Apparatus as shown:

$$E = \frac{I}{r^2}$$

E = Illumination in Lumens
 I = Intensity in cd (candela)
 r = distance in meters

PHOTOMETRY:

The Bunsen Grease Spot Photometer



Adjust to Equal Illumination

So $E_x = E_{std}$ $\frac{I_x}{r_x^2} = \frac{I_{std}}{r_{std}^2}$

Solve for I_x

Adjust the apparatus such that both sides of the screen are equally illuminated.

The illumination of the sample will equal the illumination of the standard candle which is 1 cd.

Please make your measurements and fill in the data table.

Show your method of calculations on page 2.

•

DATA TABLE

Sample Number	Unknown Light Source	r_x to screen	r_{std} to screen	$(r_x)^2$	$(r_{std})^2$	I_x
1	Normal Candle
2	Two Candles
3	Three Candles
4	Wimpy Candle
5	Gas Light

Calculations:

.

.

.

.

.

.

.

.

Critique: